

CONTENTS

- **1. BALANCE SKETCHES**
- 2. DENSITY KIT SKETCH
- **3. INTRODUCTION**
- 4. INSTALLATION

5. BALANCE OPERATIONS & FUNCTIONS

- 1) KEYBOARD DESCRIPTION.
- 2) CHANGING WEIGHING UNITS
- 3) CHANGING DATE/TIME/Sr.NO.
- 4) PRINT OPTION
- 5) RS232 INTERFACE
- 6) STORAGE OF WEIGHTS
- 7) POWER SAVING MODE.
- 8) PEAK HOLD MODE.
- 9) PIECE COUNTING MODE
- 10) SET POINT FACILITY
- 11) SIMPLE TARE MODE.
- 12) DENSITY DETERMINATION OF SOLIDS.
- 13) TEXTILE COUNTS
- 14) GSM APPLICATION
- 15) PERCENTAGE WEIGHING.
- 16) AUTOCALIBRATION WITH STANDARD WEIGHTS.
- 17) WEIGHT SLIP PRINTING.
- 6. SETUP FUNCTIONS

WARRANTY

CONTECH INSTRUMENTS LTD warrants all its products against defects in material and workmanship for a period of one year, subject to terms and conditions stated below and as further modified by warranty Amendment, in each product instruction manual. The warranty card must be registered with us within 15 days of purchase.


- 1. Warranty period will commence from the date of shipment from **CONTECH** to the original buyer.
- 2. All warranty repairs are normally done at our factory in Mumbai and our decision about faulty materials or workmanship will be final. The instrument should be sent in the original packing to our factory at the address given below. Postage /Airfreight charges both ways are to be borne by the customer.
- If any of our product is opened by any one other than our engineers or our authorised representatives, this warranty will become null and void and CONTECH will be relieved of all responsibilities as to the service and
- operation of the said product.4. This warranty will not be applicable to :
 - a) Shipping damage or damages incurred while products are in transit.
 - b) Correction of operational problems arising out of environmental conditions beyond our control.
 - c) Maintenance necessitated by customer neglect, misuse, improper operation of the instrument or equipment.
 - d) Work necessitated by damages from war, accident, fire, flood, electrical failure, vandalism or any other causalities.
 - e) Repairs due to customers failure to perform any routine maintenance prescribed in the instruction manual.
- (The routine inspection of calibration and other parameters should be done periodically by the user)
- CONTECH shall not be liable for any consequential damages nor labour loss or expense directly or indirectly arising from use of its products.,
- Amendments, assumed corollaries or statements contrary to the terms of this warranty shall not be binding on us unless they are put in writing and approved by us.
- 7. Any disputes arising out of usage of this products will be subject to Mumbai jurisdiction.
- 8. For warranty service, contact your local dealer or contact us on the below address.

Corporate Office: 301, Punit Indl. Premises, Turbhe Naka, Navi Mumbai - 400 705. Tel.: +91 22-2761 1176 / 77 / 78 / 79 / 80, 2761 8431 , 6139 3000 Fax:+91 22-2761 8421 E-mail: sales@contechindia.in / info@contechindia.in Website: www.contechindia.com

> Factory: Plot No. EL-221 TTC Indl. Area, MIDC (Electronic Zone), Mhape, Navi Mumbai-400 710. Tel.: +91 22-6194 4000 Fax: +91 22-2761 8371

Product :	Purchase Ref.:
Model No.:	Invoice No.:
Sr. No.:	Date :

Dealer's Name & Address:



INTRODUCTION

Contech[®] CAS-CAI Series weighing balances use Electro magnetic force compensation technique to measure precisely the weight of and object. The following features enable the user to suit these balances for variety of applications.



Features:

- * Multiple weighing units , Gram, Carat, Tola, , Pound , Grains, GSM, % weighing
- * Piece counting facility, up to 25 different types.
- * Storage of weights in memory and printing, up to 100 weights.
- * Power saving mode.
- * Bi-directional RS232 interface to interface with computers and printers.

1

- * Selectable baud rate.
- * Set point facility up to 2 limits.
- * Auto Power off.
- * Optional Peak Hold facility.
- * Date and time facility.
- * Multiple Print options with Sr. no., Date, Time and weight in Horizontal/Vertical mode.
- * Automatic zero tracking.
- * Density determination option.
- * Optional Battery backup facility.
- * GSM computation
- * % weighing & Calibration.
- * Weight slip printing option.

INSTALLATION

1. Unpacking:

Unpack the balance. Save the packing container for future use.

2. Electrical requirements:

The balance requires very stable power. It works on 230V AC supply with PROPER EARTHING. The power outlet used for the balance should not be shared with any other devices which draws current in inconsistent manner like Airconditioner or refrigerator etc.

3. Environmental requirements:


For best results, the balance should be placed on a level surface which is free from drafts. It should not be exposed to direct sunlight or radiated heat. The balance should not be subjected to sudden ambient temperature changes. Table used for balance should be sturdy and should not transmit vibration from other equipments and free from the movement of people. No vibration producing equipment should be operated on the same platform balance


START UP


CAS-CAI Series

Power is supplied to the scale through a 4 pin Switch mode power supply supplied along with the balance. Connect the 4 pin SMPS to the balance to a 4 pin round male connector provided at the rear panel of the balance. Insert the connector and rotate the external cover to make the connection firm and proper. Connect the 4 pin SMPS to a proper AC mains outlet with proper Earthing.

(2)







BALANCE OPERATIONS & FUNCTIONS

4. Print Option:

These balances can be attached to a serial printer for your printing needs. Print out can be programmed to suit most of the printing requirements. Note that the printer should have a serial port and baud rate of the balance and printer should be same. 2 baud rates are available are. 2400 and 4800.

key to print weight through the serial port.

Printing option and patterns are controlled by 4 SETUP parameters.

They are a) Print :	There are 4 options
I) SlogLE	 Press key to start printing weight and other details programmed as per (b), (c) and (d) below.
іі) SEABLE	- Printing is initiated when the weight kept on the pan becomes stable.
iii) ALL	 All the displayed weights are printed along with other details programmed as per (b), (c) below.
_{iv)} StorE	To print weights stored in memory along with details programmed as per (b) & (c) below.
b) Pr.Fn	Printer format. 6 printout formats are available.
iii)P-F-3 iv)P-F-4 v)P-F-5	- Serial no and weight. - Serial no, Date and weight. - Serial no., Time and weight.
	2

c) P-TYPE Print type (Horizontal or Vertical)

I) P.ESPE1 – Horizontal Details will be printed horizontally.

Sr.No.DateTimeWeight00112.05.200213:25:0023.54 g00212.05.200213:27:0523.23 g

ii) PLYPE2 - Vertical Details will be printed vertically in a slip form. For ex. Sr.No. : 001

Date : 12.05.2002 Time : 13:25:00 Weight : 24.45 g

Set the above parameters (a) to (c) to your requirements and effect printing. These parameters are available in SETUP functions.

7

BALANCE OPERATIONS & FUNCTIONS

5) BI-DIRECTIONAL RS232-INTERFACE.

Bi-directional RS-232 interface is provided in these balances to communicate with devices like computer, printer etc. The interface is provided through a nine pin D-type connector provided at the rear side of the balance. Connections are as below.

Pin 2 – RXD – Receive Data Pin 3 - TXD – Transmit Data Pin 7 – Ground.

The Serial data transmitted and received are in standard ASCII mode (+/-15V)-ASYNCHRONOUS, 8 BITS, NO PARITY, 1 STOP BIT.

Baud rate : 2400 OR 4800 SELECTABLE.

The data format for weight output is

<+/->WWWWWWWWb <bg/Ct> <CR><LF> (15 characters) where WWWWWWWW is the weight

b – blank space - 20 hex CR- Carriage Return – 0D hex LF – Line feed - 0A hex

for example, weight 85.12 g will be sent as

			_					_						1
+	0	0	0	0	8	5	1	2	b	b	g	CR	LF	

where b=black(20H), CR=carriage return (0dH) LF=line feed (0a H)

The balance could be controlled by an external device like computer with the following commands.

Z - Tares the balance.


 $\ensuremath{\mathsf{W}}\xspace\#$ - Number of times, the weight data is to be transmitted through the serial port.

can be any number from 1-9.

6. STORAGE OF WEIGHTS IN MEMORY:

Upto 100 weights can be stored in memory and recalled if required. The balance also computes the total weight of all the weights in memory. To use this option, set ACCU function to ON in SETUP functions.

a) Clear weights in memory.




Weights in memory and total weight will be cleared.

b) Storing weight in memory.


Remove all the weights from the pan and make the weight zero by pressing key.

Keep the sample on the pan and wait till the count becomes stable.



BALANCE OPERATIONS & FUNCTIONS

c) Recalling weights from memory.



(11)

(Set P	RINT option t	o" StorE	" in SETUP)
Sr.no. 1. 2. 3. 4. 5. 6. 7. 8. 9. 10.	DATE 12.10.04 12:10.04 12.10.04 12.10.04 12.10.04 12.10.04 12.10.04 12.10.04 12.10.04 12.10.04	TIME 13:30:20 13:31:05 13:31:25 13:31:50 13:32:05 13:32:40 13:32:55 13:33:20 13:33:55 13:33:55 13:34:30	WEIGHT 12.56 g 12.50 g 13.40 g 12.90 g 12.03 g 12.30 g 12.56 g 13.00 g 11.95 g 12.50 g
TOTA AVER MAXII MINIM STD. I C.V.	AGE MUM IUM	: 125.70 : 12.575 : 11.956 : 13.406 : 0.439 : 3.50%	g] g

BALANCE OPERATIONS & FUNCTIONS

7) POWER SAVING MODE

Power saving mode feature will further enhance the battery backup time by switching off the display whenever the weight displayed is zero. The balance will come out of Power saving mode when the displayed weight is not zero.


To enable this feature , refer SETUP mode.

8) PEAK HOLD MODE

Peak hold feature will enable the balance to hold the display to the maximum weight (Peak weight) displayed. When used in this mode, the balance will be continuously displaying the maximum or peak weight measured by the balance, even after the weight is removed from the pan. This feature is optional.

For using this mode, this function should be enabled in the SETUP mode.

For putting the balance into PEAK HOLD mode,




9) PIECE COUNTING MODE

Contech CAS-CAI series balances can be used for piece counting purposes also. Piece calibration of 25 items can be stored in memory. Accuracy of piece counting depends on the uniformity in weight of the items and the sample size used for piece calibration. Better the weight uniformity and more the sample size, better will be the accuracy.

Use SETUP function to select proper piece counting mode before using. There is an option to select 1, 10, 25 or none piece counting modes. See relevant section in SETUP functions for more details.

Selection of piece counting memory (item)



Piece Calibration

Select proper piece counting memory as mentioned above.

Make the Weight read zero by pressing the key. Weight of any

(13)

 $\widehat{}$

key.

 $container/bags \ used should also be made zero by pressing$

Keep known number of pieces on the pan.

Press (F) key, the balance will display


PC xxxx

BALANCE OPERATIONS & FUNCTIONS

Use equal and keys to

keys to change the number of pieces and make it

to the pieces kept on the pan.



1

Balance does piece calibration and displays the number of pieces on the pan, say $% \left({{{\mathbf{x}}_{i}}} \right)$

250

10) OPERATING IN SIMPLE TARE MODE:

In SIMPLE mode, balance can be operated in its simplest mode. In this mode,



to make the weight zero.

All other functions will be disabled

(14)

11) SET POINT FACILITY

This facility enables the user to set up to 2 weights for comparison with the current weights to activate different events. This feature is controlled by 3 SETUP functions.

1. F{LL - Make it El -no 2. SEF-BF - Make it SEL-7 for 1 or 2 Set Points 3. F_EYPE swapping - Make it F-TYPE 1 or F-TYPE 2 for event depending on the requirement. For setting weights, Press (SET UP key, and immediately followed by key. the balance will show FILL- I and followed by XXXX.XX Where XXXX.XX is the current set limit - 1. Use (A) and (keys to change the limit. To discard changes, Press () & keys simultaneously or (ENTER) Press key to save. The balance will show FILL-2 (If fill option is on and limit is set to SET-2) XXXX.XX For some time and will display Where XXXX.XX is the current set limit. Use (RINT) and (V) keys to change the limit. To discard changes, Press (& & kevs simultaneously or Press key to save. (15)

BALANCE OPERATIONS & FUNCTIONS

13) DENSITY DETERMINATION OF SOLIDS

This is an optional feature, valid only if this facility is available in the product supplied. Balance calculates the density of solids based on Archemedis principle and displays it. The sample is weighed in both Air & Water and the balance calculates the density. There is a provision for entering the water temperature, so that necessary correction is applied for calculating density.

There are 2 types of density measurements.

1. Normal mode

D = Weight in Air D = x DL x DW (Weight in Air - Weight in Liquid)

D – DensityDL – liquid density (Programmable)DW – Density of Water at measuring temperature.(If water is used)

For DW, Water temperature is programmable from 10 to 50 deg.C

DL is also programmable.


If DL is set to 1.0000 then water temperature correction will be applied.

If DL is set to a value other than 1.000 then water temperature correction will not be applied.

For example Weight in Air = 8.323g, Weight in water = 6.222g Water temperature = 25deg. Liquid density = 1.0000

Density = 3.9497 g/cc

(16)



13) TEXTILE COUNTS

CAS/CAI Series balances directly calculate and display textile counts like cotton count, denier count, Tex Count, Worsted Count and sliver. This is an optional feature. The following details are valid, only if this feature is available in the balance. Enable this option in SETUP functions. (See Setup functions for more details). Length of the sample for the above counts are as follows.

Cotton count - 120 yd. dEntEr Denier count - 90 m. EEX Tex count - 10 m Gor SEEd Worsted count - 10 yd. SLiger Sliver - 10 yd.	
Use (RINT) and (V) keys to change the unit ,	
For eg. Cotton	
Balance will display cotton count directly. It is also possible to analyse the samples for consistency using statistical method.	
For using statistical method , enable the SERE and RECU options in SEE UP	
functions. Also set the PRINT option to StorE	
For beginning the tests, clear the memory by the following.	
Press key, the balance will display	
Press (Ref) key immediately, balance will display StorE	
Press 🛞 key, balance will display	
Press key to clear the memory.	
Store counts of all the samples to be tested one by one as per the following. Remove all the weights from the pan and make the weight zero by	
pressing key.	
(19)	

BALANCE OPERATIONS & FUNCTIONS

Keep the sample on the pan and wait till the count becomes stable key, the balance will display Press PRINT key immediately, balance will display Press StorE key to store the count in memory. Press Store all the counts in memory with above procedure. For viewing sample data and other statistical details. do the following. Press key, the balance will display StorE Press key immediately, balance will display Press key, balance will display XX (displays the sampleno.) (🔻 $\Box\Box$ XXX.XX (displays the count) followed by and (🔻 keys to see all sample data. After the balance displays Press the last data, it will display Total count and total number of samples. Then it displays AU xxx.xx (Average count) Press key, balance displays (Maximum count) XXX.XX T ΠĂ (Minimum count) Press key, balance displays XXX.XX . _ $\cap \cap$ Press key, balance displays XXX.XX (Standard deviation) V 58 key, balance displays xxx.xx (Coefficient of Variation) Press T 68 While viewing the above, Press PRINT key to print all the results on a serial printer. Sample printout will appear as below. COUNT Sr.no. 12.56 CC 1. 12.50 CC 2. 3. 13.40 CC TOTAL : 125.70 12.90 CC 4. **AVERAGE : 12.57** 5. 12.03 CC MAXIMUM : 11.95 12.30 CC 6. MINIMUM : 13.40 7. 12.56 CC 8. 13.00 CC STD. DEV. : 0.439 9. 11.95 CC C.V. : 3.50% 10. 12.50 CC (20)

14) GSM APPLICATION

CAS-CAI series balances can be used to determine GSM(Grams per Sq.meter). The following description is valid only if this feature is available in the balance and enabled. The balance displays directly the GSM value of fabric or paper of specified area. Enable GSM feature in the SETUP functions. 5 pre-programmed area (rectangular or round) are available in the balance. Standard areas are:

1. 5 X 5 cm 2 10 X 10 cm 3. 20 X 20 cm 4. 25 X 20 cm 5. 25 X 25 cm

Select GSM mode:

PRINT V keys to change till the balance displays -05and Use

Balance now enters GSM mode.

For changing the area, Press key, the balance displays

Press key immediately, balance displays XX.XX

Repeat the procedure till desired area is selected.

15) PERCENTAGE WEIGHING

Selection of percentage weighing mode.

Press keys till the balance displays or PRINT

PERCENTAGE weighing function is used to determine % weight gain/loss. Any weight within the capacity can be set to be 100%. The balance displays the weight gain/loss in % of the original weight. This function must be enabled using SETUP function. This feature is very useful in determining % loss/gain in moisture in food/tea/pharmaceutical industries.

PFee

For eg. Keep a 25g weight on the pan.

F. key, the balance calibrates 25g as 100% and Press

displays 100.0

After this any change in weight on the pan will be indicated as % of original weight(25g), till the balance is calibrated for100% with another weight.

(21)

BALANCE OPERATIONS & FUNCTIONS

16) AUTOCALIBRATION WITH STANDARD WEIGHTS

CBB-Series balances can be calibrated for weight with standard mass. Balances can be calibrated with 100g, 200g, 500g and 1000g weights depending on the models.

Models upto 1kg can be calibrated with 100g,200g,500g, 1000g weights Models upto 600g can be calibrated with 100g,200g,500g weights Models upto 350g can be calibrated with 100g,200g weights Models upto 300g can be calibrated with 100g,200g weights Models upto 220g can be calibrated with 100g,200g weights Models upto 125g can be calibrated with 50g, 100g weight.

Enable calibration function in SETUP mode before attempting to calibrate the balance.(Refer SETUP functions for more details). This function should not be made available to the end user, if there is any restriction in usage of this function.

CALIBRATING THE BALANCE

Use only good calibrated weights for performing auto calibration. Press 🕥 key to make the weight read zero.

Keep the standard mass on the pan and wait for it to become stable.




Balance will display And will subsequently display

For setting balance calibration back to default factory setting,


key and immediately press (F) key. Press







SETUP FUNCTIONS



SETUP FUNCTIONS

SETUP functions control the basic operation of the balance. There are 29 parameters, which can be set by the user to suit the requirements. The following are the parameters.

MENU NAME ¹ FRCL	FUNCTION Factory setting	OPTIONS	DESCRIPTION To select factory set parameters.
2. PA-tS	Piece counting modes	PR-t-1 PR-t-10 PR-t-25 00 PR-t	To select single piece counting memory. To select 10 piece counting memory. To select 25 piece counting memory To disable piece counting.
3. 6808	Baud rate setting	695400 694800	To select 2400 baud rate To select 4800 baud rate
4. Print	Print modes set.	SingLE	Send stable weight through serial port
			when key is pressed.
		ՏեԵ	Send weight through serial port
		ALL StorE	Every time balance reading becomes stable. Send weight continuously. Send stored weights through serial port.
(Ability	Auto zero setting of the balance to remain when there is no weight ban.	0-8 8-1 8-2 8-3	Autozero disabled. Autozero to half accuracy of balance. Autozero to full accuracy of balance. Autozero to twice the accuracy of balance.
6. ACCU	Weight storage mode	Ac-no Ac-YES	Weight storage disabled Weight storage enabled
7. FILL	Fill mode option	FL-oFF FL-on	Fill mode disabled Fill mode enabled
		25	

SETUP FUNCTIONS

^{8.} -[AL-	Auto calibration	(AL-oFF (AL-on	Autocalibration disabled Autocalibration enabled
9. Pound	Pound weighing	Pnd-off Pnd-on	Pound weighing disabled Pound weighing enabled
^{10.} (A-AL	Carat weighing	(rt-off (rt-on	Carat weighing disabled Carat weighing enabled
^{11.} Łola	tola weighing	tol-off tol-on	Tola weighing disabled Tola weighing enabled
12. GrAin	Grain weighing	GrAin-off GrAin-on	Grain mode disabled Grain mode enabled
13. PErC	% weighing	PEr[-off PEr[-on	% mode disabled % mode enabled
^{14.} 656	GSM	657-0FF 697-09	GSM mode disabled GSM mode enabled
15. PS 7	Power saving mode	PSā-off PSā-on	Power saving mode disabled Power saving mode enabled
16. P-HoLd	Peak hold mode	РНЕ-0 ⁶⁶ РНЕ-00	Peak hold mode disabled Peak hold mode enabled
		(26)	

SETUP FUNCTIONS

17. 33-5EE Third decimal mode	РЭд -oFF РЭд -on	Third decimal mode disabled Third decimal mode enabled
18. ELE Title printing	ELE -off ELE -on	Title printing disabled Title printing enabled
19. P-LYPE Select printing mode	P-EYPE (P-EYPE 2	Horizontal Printing mode Vertical Printing mode
20. Pr.,FmL Select print formats	Pr.F { Pr.F 2 Pr.F 3 Pr.F 4 Pr.F 5 Pr.F 6	Only weight Sr.no. Weight Sr.no., Date, Weight Sr.no., Time, Weight Sr.no., Date, Time, Weight Weight slip Printing
21. SEE-PE Select set point mode	SEE-1 SEE-2	Single set point 2 set points.
22. F-LYPE Select set point type	F-EABE (E-FABE S	+ve Logic outputs for fill application -ve Logic outputs for fill application
23. dEn ^S Density mode	dEnS-oFF dEnS-on	Density mode off Density mode on.
24. Count Textile count mode.	Count-off Count-on	Textile Count mode off Textile Count mode on.
25. SERE Statistical report mode	SEAE-oFF SEAE-on	Statistical report off Statistical report on.
	27	

SETUP FUNCTIONS

PArts	It parameters.	To coloct single piece counting memor
<u>nnrco</u> bRUd		To select single piece counting memor
	6d4800	To select 4800 baud rate
Print	Single	Send stable weight through serial port
RUEo-O	R-1	Autozero to half accuracy of balance
REEU	Ac-no	Weight storage disabled
FILL	FL-off	Fill mode disabled
(AL	(RL-off	Autocalibration disabled
Pound	Pnd-off	Pound weighing disabled
(ArAt	Crb-off	Carat weighing disabled
tolЯ	tol-off	Tola weighing disabled
grflin	GrBin -off	Grain mode disabled
%	PEr(-off	% mode disabled
656	657-off	GSM disabled
P57	PSA-off	Power saving mode disabled
P-Hold	рна-орр	Peak hold mode disabled
3d-58E	639 -0 ₂₅	Third decimal mode disabled
FIFTE	LLE -off	Title printing disabled
P-EYPE	P-E9PE (Horizontal Printing mode
Pr,Filt	Pr.,F (Only weight
SEE UP	SEE-1	Single set point
F-EYPE	F-ŁYPE (+ve Logic outputs for fill application
Count	(ount-off	Textile counts off
d8~£9PE	dEnS-oFF	Density type 1
SEAE	SERE-OFF	Statistical reports off-