

BALANCE ACCESSORIES

BATTERY BACKUP

Battery backup unit is an online uninterrupted power supply providing adequate voltages for the operation of the balance in case of power failure. Battery Backup unit incorporates sealed maintenance free lead acid batteries with a charger and under voltage and over voltage protection circuit.



CONTECH[®]
Instruments Ltd.
AN ISO 9001, 14001 & 45001 COMPANY

Where innovation never stops.....

301, Punit Indl. Premises, Turbhe Naka, Navi Mumbai - 400 705.

Tel. : +91 22 6139 3000 (100 lines)

E-mail : sales@contechindia.in

Website: www.contechindia.com

FACTORY / NABL LAB



Plot No. EL-221 TTC Indl. Area, MIDC (Electronic Zone),
Mhape, Navi Mumbai - 400 710.