

BALANCE ACCESSORIES

SWITCH MODE POWER SUPPLY (SMPS)

Contech scales are supplied with SMPS which provides required DC voltages for the operation of the scale circuit. SMPS operates on a high frequency and hence the circuit becomes efficient & compact. The input AC power voltage span is 100-270V, AC, 50/60Hz and hence with SMPS, the scale can operate on both 110V and 220V AC supply.



CONTECH[®]
Instruments Ltd.
AN ISO 9001, 14001 & 45001 COMPANY

Where innovation never stops.....

301, Punit Indl. Premises, Turbhe Naka, Navi Mumbai - 400 705.

Tel. : +91 22 6139 3000 (100 lines)

E-mail : sales@contechindia.in

Website: www.contechindia.com

FACTORY / NABL LAB



Plot No. EL-221 TTC Indl. Area, MIDC (Electronic Zone),
Mhape, Navi Mumbai - 400 710.